

Waste Not, Want Not: Optimising Wastewater Surveillance of Human Enteroviruses.

Ms. Teesta Maulik¹

Dr. Ben Bellekom¹

Dr. Alexander Shaw¹

¹School of Public Health, Imperial College London

Contact: teesta.maulik22@imperial.ac.uk

01 Introduction

- Enteroviruses, one of the most prevalent genera of viruses globally, can pose a significant threat, causing up to 80% of viral meningitis cases worldwide.¹ Early detection of dangerous serotypes in wastewater is crucial for implementing control measures.
- They are RNA viruses with an outer protein coat made of four polypeptides including VP1, VP2, VP3 and VP4.²
- The VP1 region is typically used to identify the serotype of the virus.²
- Polymerase chain reaction (PCR) assays can be used to amplify target regions in the genome for identification via sequencing. The PCR specificity relies on short DNA fragments known as primers.
- Two common primer sets (NCR_CRE and MM_EV) used to amplify a large part of the enterovirus genome (which includes the VP1 region) have performed poorly in several wastewater surveillance labs worldwide.³
- This project aims to compare the efficacy of alternative primers in detecting enteroviruses For1_5NCR/Rev7_5NCR and For2_5NCR/Rev6_5NCR (designed by Dr. Alexander Shaw) against current primers using identical samples. These primers do not target the VP1 region and hence, would only identify the species of enteroviruses.




Figure 1: Primer positions on the Sabin AY184219 Reference (approximation).

Red = For1_5NCR and Rev7_5NCR

Purple = For2_5NCR and Rev6_5NCR

Blue = 5'-NTR / MM_EV and CRE / MM_EV

02 Method

Viral RNA Extraction from wastewater using the MagMax Viral RNA isolation kit

Nested Reverse-Transcription Polymerase Chain Reaction (RT-PCR). Then, verification using Gel Electrophoresis and TapeStation®

Amplicon clean-up using the AMPure XP beads. End-prep & dA-tailing, then ligation of barcodes using the Native Barcoding kit 96 (SQK-NBD114.96)

Preparation of, and sequencing amplicons on the ONT® MinION Flow Cell. Data was read-mapped via minimap2 and analysed on Geneious Prime®.

03 Results

RT-PCR was performed on positive/negative controls and wastewater replicates. Amplicon sizes were verified via gel electrophoresis (see Figure 2 and 4) and TapeStation® (see Figure 3). After which, the PCR products were sequenced. An arbitrary threshold of over 50 reads was interpreted as positive detection. The resulting data (see Table 1 and 2) shows the average number of reads of each enteroviral species detected per primer set for wastewater amplicon replicates and the negative control.

1. 5'-NTR-CRE primers: No species detected.
2. MM_EV primers: Detected species A and B.
3. For1_5NCR/Rev7_5NCR, For2_5NCR/Rev6_5NCR, and nested primers: High read counts for species A, B, and C.
4. Negative control: Positive for species C with new primers.




Figure 2: Agarose gel electrophoresis of PCR products with 1% agarose gel with 1kb ladder. Lanes show amplification products from wastewater samples using NCR_CRE, MM_EV, and For2_5NCR/Rev6_5NCR primers sets. The NCR_CRE and MM_EV also show amplification products with positive and negative controls.




Figure 3: TapeStation analysis of PCR products with 1kb ladder. Lanes show amplification products from wastewater samples using For1_5NCR/Rev7_5NCR primers sets. The NCR_CRE and MM_EV also show amplification products with positive and negative controls.




Figure 4: Agarose gel electrophoresis of PCR products with 1% agarose gel with 100bp ladder. Nested RT-PCR with For2_5NCR/Rev6_5NCR primers on For1_5NCR/Rev7_5NCR amplicons.

Primer Sets / Species Type	5'-NTR_CRE	MM_EV	For1_5NCR and Rev7_5NCR	For2_5NCR and Rev6_5NCR	For2_5NCR/Rev6_5NCR Nested into For1_5NCR/Rev7_5NCR
A	8	155	690	560	10866
B	1	190	1074	1018	2217
C	16	9	8472	42939	12884
D	0	0	1	8	26

Table 1: Following PCR product verification, the table shows the average number of reads of each enteroviral species detected per primer set for wastewater amplicon replicates.

Primer Sets / Species Type	5'-NTR_CRE	MM_EV	For1_5NCR and Rev7_5NCR	For2_5NCR and Rev6_5NCR	For2_5NCR/Rev6_5NCR Nested into For1_5NCR/Rev7_5NCR
A	0	0	0	1	2
B	0	0	1	0	0
C	0	9	700	37684	28125
D	0	0	0	0	0

Table 2: Following PCR product verification, the table shows the average number of reads of each enteroviral species detected per primer set for the negative control.

04 Discussion

- The Majumdar et al. study utilised 5'-NTR-CRE and MM_EV primer sets for enterovirus species identification.³ However, this project revealed non-specific binding and failed detection in positive controls, corroborating recent observations from other laboratories. This discrepancy may be attributed to alterations in reagent composition.
- Figure 1 illustrates that 5'-NTR-CRE amplifies a ~3.8kb genomic region, while the new primers target a ~600bp segment. Longer amplicons necessitate intact viral RNA for RT-PCR, but RNA's inherent instability complicates accurate amplification of extended sequences. Additionally, longer target regions increase the likelihood of DNA depurination and damage during PCR.⁴
- Observed differences in read counts may result from incomplete normalisation or variations in MinION sequencing depth, rather than solely differential PCR efficiency.
- As expected, nested PCR yielded the highest sequencing reads due to its two-step amplification process. The short target region and nested approach enhanced sensitivity and specificity, resulting in increased target sequence yield and consequently higher read counts.
- Negative controls across all experiments exhibited contamination with the positive control species, suggesting potential carry-over.
- The use of Coxsackievirus A20 as a positive control has implications for method validity. The absence of species D enterovirus detection warrants further investigation to determine whether it reflects a true absence or insufficient primer sensitivity for this species.

05 Conclusion and Next Steps

- This experiment shows that the current primers used for gold-standard diagnostics are less sensitive and alternative primers provide a method to detect enteroviral species diversity in sewage.
- Further work must be done on isolates to show more accurate results and method validity.
- Consider further work to explain why the NTR-CRE and MM-EV primers may not be working.
- The alternative primers bind outside the VP1 region, allowing only species detection. Future work should develop VP1-specific primers for serotype identification.

06 Acknowledgements and References

Sincere thanks to the Laidlaw Foundation for generously sponsoring this project, without which I would not have been able to pursue this opportunity. I would like to express my heartfelt thanks to my project supervisor, Dr. Alexander Shaw for his consistent guidance and patience. Thank you to the Vaccine Epidemiology Research Group (VERG) and Professor Nicholas Grassly at Imperial College London for hosting me. Thank you to Dr. Ben Bellekom, Ms. Catherine Troman, Ms. Shannon Fitz and Dr. Jaspreet Mahindroo for all your support through this project.

1. Greenlee, J. (2014) 'Nervous System Complications of Systemic Viral Infections', in Aminoff's Neurology and General Medicine (Fifth Ed.) pp. 858
2. Swiss-Prot group of the SIB Swiss Institute of Bioinformatics. Available at: <https://viralzone.expasy.org/97> Accessed: 06 August 2024
3. Majumdar, M. and Martin, J. (2018). Detection by Direct Next Generation Sequencing Analysis of Emerging Enterovirus D68 and C109 Strains in an Environmental Sample From Scotland. Frontiers in Microbiology, 9. doi:https://doi.org/10.3389/fmicb.2018.01956.
4. Qiagen. (no date). PCR amplification: Long PCR products. Available at: <https://www.qiagen.com/us/knowledge-and-support/knowledge-hub/bench-guide/pcr/introduction/amplification-of-long-pcr-products>. Accessed: 10/10/24