


Open Sesame: Can We Make a Protein Cage Open at the Right Time?

Engineering a pH-Responsive Protein Nanocage for Controlled Cargo Release

Kamalini Rajesh | Supervisors: Dr. Stefanie Frank and Dr. Yuqian Ou

1. Background


- Encapsulins are **self-assembling bacterial protein nanocages** that can protect molecular cargo, but their **high stability can limit controlled release** [1].
- The **mildly acidic tumour microenvironment** provides a **potential trigger** for pH-responsive cargo release [7].
- Here, we engineered the **pH-sensitive peptide GALA** into the **E-loop of T=4 *Quasibacillus thermotolerans* encapsulin (QtEnc)**.
- Under **acidic conditions**, GALA changes conformation and is hypothesised to disrupt subunit interactions, promoting cage disassembly.




2. Aims

- Express and purify wild-type EncStrep and T4GALA EncStrep in *E. coli*
- Characterise the pH-dependent assembly behaviour of T4GALA
- Generate T4GALA EncStrep co-expressing GFP-CLP as a model cargo system
- Assess GFP cargo loading and investigate evidence of pH-triggered release

3. Engineering the pH responsive system




4. Overview of Methods




5. Results


5.1 Size Exclusion Chromatography (SEC)




5.2 SDS-PAGE Analysis



5.3 Native PAGE: pH Dependent Assembly



5.4 T4GALA EncStrep + GFP-CLP: Cargo Loading



6. Conclusions

- T4GALA retained assembly under neutral and basic conditions, with reduced cage-associated signal under more acidic conditions.
- T4GALA EncStrep + GFP-CLP was successfully generated, and GFP-CLP co-purified with the encapsulin system.
- GFP fluorescence was associated with cage species, supporting GFP cargo association.
- However, clear pH-triggered release of GFP from the cage was not observed and could not be conclusively demonstrated.
- An unexpected lower-molecular-weight species was observed with T4GALA and requires further characterisation.

7. Future Work

- Identify the lower-MW species (e.g. mass spectrometry)
- Quantify encapsulated vs released GFP (e.g. fluorescence spectroscopy)
- Test Tumour-relevant pH range (6.5-7.4) using a controlled buffer system
- Optimise expression and purification conditions to reduce subunit heterogeneity

Scan to access the full report and references.



Funding from the UK Engineering and Physical Sciences Research Council (EPSRC) for the Manufacturing Research Hub for a Sustainable Future (VaxHub Sustainable), co-directed by University College London and the University of Oxford with UK university partners is gratefully acknowledged (Grant Reference: EP/X038181/1)